


MLRB_LHTWIN_SS

Bariatric Drop Down Grab Rail (Twin), Left Hand Nurse Call, Right Hand Toilet Roll Holder, Drop Down Leg



SPECIFICATION

Material	#304 Grade Stainless Steel 32mm with solid tube reinforcement
Finish	Satin Stainless Steel
Tubing	Solid Tube Reinforcement
Dimensions	830mm from the wall x 270mmH
Mounting Bracket & Brace	#304 Grade Stainless Steel, 10mm
Mounting Bracket	270mmH x 78mmW x 10mmThick Plate
Mounting Holes	6 x 5.5mmØ
Drop Down Leg	600mm (The Drop Down Leg has approx 10mm of extension to allow for uneven surfaces)
Nurse Call Housing	Fixed LH/RH Options
Toilet Roll Holder	Toilet Paper Holder changes easily with a spring loaded mechanism
Features	Unit is Seamless Welded & Polished Click-lock system when the grab rail is in the upright position

NOTES

- Special Order Product - Please contact Metlam Australia for Details**
- For Metlam Australia Care and Maintenance Guidelines, Installation Instructions, Product Specifications, and Warranty Details, please refer to the Product User Guide Tab.

 For more in depth information, visit the product online:
https://metlam.com.au/category/grab-rails/MLRB_LHTWIN_SS