


MLRB_RHTWIN_SS

Bariatric Drop Down Grab Rail (Twin), Right Hand Nurse Call, Left Hand Toilet Roll Holder, Drop Down Leg



SPECIFICATION

Material	#304 Grade Stainless Steel 32mm with solid tube reinforcement
Finish	Satin Stainless Steel
Tubing	Solid Tube Reinforcement
Dimensions	830mm from the wall x 270mmH
Mounting Bracket & Brace	#304 Grade Stainless Steel, 10mm
Mounting Bracket	270mmH x 78mmW x 10mmThick Plate
Mounting Holes	6 x 5.5mmØ
Drop Down Leg	600mm (The Drop Down Leg has approx 10mm of extension to allow for uneven surfaces)
Nurse Call Housing	Fixed LH/RH Options
Toilet Roll Holder	Toilet Paper Holder changes easily with a spring-loaded mechanism
Features	Unit is Seamless Welded & Polished Click-lock system when the grab rail is in the upright position

NOTES

- Special Order Product - Please contact Metlam Australia for Details**
- For Metlam Australia Care and Maintenance Guidelines, Installation Instructions, Product Specifications, and Warranty Details, please refer to the Product User Guide Tab.

[For more in depth information, visit the product online:
https://metlam.com.au/category/grab-rails/MLRB_RHTWIN_SS](https://metlam.com.au/category/grab-rails/MLRB_RHTWIN_SS)